EN Increasing the predictability of ceramic veneers by using a resin cement with an advanced polymerization system (APS)

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32-year-old male patient.

COMPLAINT

Dissatisfaction with his smile due to several interdental spaces in the upper anterior teeth.

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EN CLINICAL/RADIOGRAPHIC EXAM

The presence of diastemata was noted in the region of the central, lateral and upper canines.

TREATMENT

Due to the search for an immediate solution to the patient's case and a good oral hygiene condition, ceramic veneers were chosen to improve the smile with: four veneers (12 to 22) and two ceramic fragments (13 and 23).

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2a and 2b. Shade taking.
2a y 2b. Tomada de color.

3a and 3b. Mock-up installation with bisacrylic resin according to waxing done in a model. This step allows the patient to preview the final aesthetic result, with the possibility of minor adjustments and corrections. With the approval of the new smile aesthetics by the professional and the patient, the teeth preparation phase began.

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4a and 4b. Dental preparation guide in silicone created on the waxed model and adapted to the teeth already prepared, evidencing the amount of wear obtained. In this case, the modified prepless technique was used, where a cervical preparation of approximately 0.2mm was made on the gingival line with a round-tipped cylindrical diamond tip (KG Sorensen) mounted on a contra-angle multiplier with a rotation of around 50,000rpm, allowing precision and delicacy. In the proximal ones, a preparation of 0.2mm was also carried out with a long conical diamond tip. The preparations were carried out on the mock-up in position with the aid of "silicone guides" that allowed the precise and conservative wear necessary for the ceramics. These preparation guides allow checking of the worn thickness in the different regions, avoiding the unnecessary removal of healthy structure.

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
5a to 5c. Insertion of the retractor cord continuously and mold obtained in addition silicone (double-molding technique), with a true copy of the cervical finish lines of the preparations.

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6a and 6b. No-prep veneers and ceramic fragments on the plaster model, in vestibular and palatal view.

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
7. Dental prophylaxis to promote an adequate substrate for the adhesive procedure.

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8. Dry proof of pieces for an adaptation check.

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
9. Allcem Veneer APS system with Try-in cements and pastes and Ambar APS adhesive.

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10. Test with Try-in pastes of different colors: OW (opaque white) on tooth 21 and A3 on tooth 11. The influence of cement color on the final aesthetic result of the tooth is evident, so the professional should have an adequate range of colors to choose from at their disposal. In this case, after the consent of the patient, the color TRANS (translucent) was chosen for cementation.

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11a and 11b. Etching of the pieces with hydrofluoric acid at 10% (Condac Porcelana) for 20 seconds. It could also be etched with Condac Porcelana 5%. After careful washing with water, application of 37% phosphoric acid (Condac 37) for 30 seconds for removal of the silica deposits created during etching, drying and application of silane (Prosil) for 60 seconds, one must wait for 1 minute for the solvent to evaporate.

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12a and 12b. At the cementation stage, teeth receive a continuous retraction cord to expose the cervical finish lines properly, control the moisture from gingival sulcus and facilitate the removal of excess cement within the gingival sulcus after light curing.

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13a to 13c. Dental preparations being etched with 37% phosphoric acid (Condac 37) for 15 seconds. After washing for 30 seconds and drying, an opaque and frosted appearance of the enamel is obtained.

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14a and 14b. Application of dental adhesive (Ambar APS), with a practically colorless appearance that does not interfere in the final color of the work.

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